







| SCHEDA TECNICA | TECHNICAL FILE | |
|---|---|---|
| GAMMA DI PRODOTTO: Sanitari ed accessori | PRODUCT RANGE: Sanitaryware and accessories | CODICE /CODE CL-TB20-W |
| DESCRIZIONE: Vaso sanitario filoparete "Giada" con altezza maggiorata | DESCRIPTION: 'Giada' floor standing WC with increased height |  |
| CARATTERISTICHE DEL PRODOTTO: <ul style="list-style-type: none"> Vaso sanitario filoparete in ceramica bianca con altezza maggiorata (H50 cm, per aiutare chi ha difficoltà a sedersi e ad alzarsi). Con scarico universale. Disponibile sedile con coperchio (codici: CL-TSL7-W e CL-TSL8-W, versione slim). | PRODUCT FEATURES: <ul style="list-style-type: none"> Floor standing WC in white ceramic with increased height to help those who have difficulty sitting down and getting up. With universal drain. Available with WC seat cover (codes: CL-TSL7-W and CL-TSL8-W, slim version). |  |
| MATERIALI: Ceramica | MATERIALS: Ceramic | DIMENSIONI /DIMENSIONS |
| FINITURA: - | FINISH: - |  |
| PRESTAZIONI MECCANICHE: - | MECHANICAL PERFORMANCE: - |  |
| CERTIFICAZIONI: Dichiarazione di conformità CE alla direttiva Europea dei dispositivi medici. | CERTIFICATIONS: CE Declaration of conformity of the European Directive for Medical Device 93/42 CEE. |  |
| NORME APPLICABILI: UNI EN 12182 UNI CEI EN ISO 14971 UNI EN 1041 | APPLICABLE NORMS: UNI EN 12182 UNI CEI EN ISO 14971 UNI EN 1041 |  |
| MANTENIMENTO: Controllare periodicamente il corretto fissaggio. Effettuare un adeguata pulizia con acqua e sapone: non utilizzare prodotti abrasivi. | MAINTENANCE: Periodically check that fittings are conform to initial mounting. Clean with water and non-aggressive detergent. Do not use abrasive or corrosive products. | |