



Ⓜ **Istruzioni di lavaggio:** Detergenti neutri e non abrasivi - Autoclave a 85 °C per un minuto - 
 Ⓔ **Washing instructions:** Detergent neutral and non-abrasive - autoclave at 85 °C for one minute - 
 Ⓡ **Instructions de lavage:** Détergent neutre et non abrasif - autoclave à 85 °C pendant une minute - 
 Ⓔ **Instrucciones de lavado:** Detergente neutro y no abrasivo - autoclave a 85 °C durante un minuto - 
 Ⓓ **Waschanleitung:** Waschmittel neutral und nicht-abrasive - Autoklaven bei 85 °C für eine Minute

